(90 × 210mm Size)

FUJIFILM

Wako

〈For Research Use Only〉 Code No. 294-63601 (750 Tests)

LabAssay™ NEFA (ACS-ACOD method)

〔Introduction〕

Non-esterified fatty acid (NEFA) in serum binding albumin, is used as an important energy source of peripheral tissues.

LabAssay™ NEFA is the reagent kit for assay of NEFA based on an enzymatic method using 3-methyl-*N*-ethyl-*N*-(β-hydroxyethyl)-aniline (MEHA) as a violet color agent. The kit is used for the quantitative determination of NEFA in mouse serum. It is a simultaneous multi-sample assay format using a microplate, but measurements can also be made using a test tube.

〔Kit contents〕

(1)	Chromogen Reagent A After reconstitution, Coenzyme A (CoA) 0.73mmol/L ATP disodium salt 4.5mmol/L Acyl-CoA synthetase (ACS) 0.27units/mL 4-Aminoantipyrine 1.5mmol/L Ascorbate oxidase (AOD) 2.7units/mL	For 10mL	6vials
(2)	Solvent A (for Chromogen Reagent A) (Phosphate buffer, pH 7.0)	65mL	1vial
(3)	Chromogen Reagent B After reconstitution, Acyl-CoA oxidase (ACOD) 5.5units/mL Peroxidase (POD) 6.8units/mL	For 20mL	6vials
(4)	Solvent B (for Chromogen Reagent B) (3-methyl- <i>N</i> -ethyl- <i>N</i> -(β-hydroxyethyl)- -aniline (MEHA) 1.2mmol/L	130mL	1vial
(5)	NEFA Standard Solution (Oleic acid 1mEq/L)	10mL	1vial

〔Materials and apparatuses to be prepared〕

- 96wells microplate (transparent type)
- Micropipette
- Incubator maintaining at 37°C *
- Plate mixer*
- Microplate reader with 550nm wavelength filter
(* if the microplate reader is not equipped.)

〔For test tube method〕


- Test tube
- Pipette
- Spectrophotometer or colorimeter with 550nm wavelength filter

—1/10—


〔Assay principle〕

Non-esterified fatty acid (NEFA) in the sample is converted to Acyl-CoA, AMP and Pyrophosphoric acid (PPi) by the action of Acyl-CoA synthetase (ACS), under coexistence with Coenzyme A (CoA) and Adenosine 5'-triphosphate disodium salt (ATP). The obtained Acyl-CoA is oxidized and yields 2,3-*trans*-enoyl-CoA and hydrogen peroxide by the action of Acyl-CoA oxidase (ACOD). In the presence of peroxidase (POD), the hydrogen peroxide yields a blue purple pigment by quantitative oxidation condensation with 3-methyl-*N*-ethyl-*N*-(β-hydroxyethyl)-aniline (MEHA) and 4-aminoantipyrine. NEFA concentration is obtained by measuring absorbance of the blue purple pigment.

ACS

ACOD

POD

〔Preparation of reagents to be used〕

① Color Reagent A :

Prepare Reagent A by mixing 1 vial of each Chromogen reagent A (for 10mL) and Solvent A (10mL). After reconstitution, the solution should be stored below 25°C and used the same day, or stored at 2-10°C and used within 5 days.

② Color Reagent B :

Prepare Reagent B by mixing 1 vial of each Chromogen reagent B (for 20mL) and Solvent B (20mL). After reconstitution, the solution should be stored below 25°C and used the same day, or stored at 2-10°C and used within 5 days.

③ Standard solution (Microplate method)

Standard solution is prepared by dilution of the provided Standard.

No.	NEFA Standard	Distilled or deionized water	Sample volume	NEFA Concentration
1	1.0mL	1.0mL	4 μL	0.50 mEq/L
2	Undiluted	—	4 μL	1.00 mEq/L
3	Undiluted	—	8 μL	1.97 mEq/L*1

*1 The test sample volume is usually 4 μL, but 8 μL is taken in this case. The NEFA concentration must be corrected accordingly as indicated in the table above.

—2/10—

[Procedure]

(1) Assay in a microplate

Perform the assay in the wells according to the following table scheme.

	Test	Standard	Reagent Blank	Sample Blank* ²
Color Reagent A	80 μL	80 μL	80 μL	80 μL
Specimen	Serum 4 μL	Standard 4 μL	Distilled or Deionized water 4 μL	—
Mix well and incubate at 37°C for 10 min.				
Color Reagent B	160 μL	160 μL	160 μL	160 μL
Specimen	—	—	—	Serum 4 μL
Mix well and incubate at 37°C for 10min.				
After cooling the reaction solution to room temperature, measure the absorbance at 550 nm of the test sample and standard solution with the blank solution as the control within the next 30min.				

*2 Sample blank is not required for common samples, but required for chyle or hemolyzed sample.


(2) Assay in a test tube

Perform the assay in a test tube according to the following table scheme.

	Test	Standard	Reagent Blank	Sample Blank* ²
Color Reagent A	1.0 mL	1.0 mL	1.0 mL	1.0 mL
Specimen	Serum 0.05 mL	Standard 0.05 mL	Distilled or Deionized water 0.05 mL	—
Mix well and incubate at 37°C for 10 min.				
Color Reagent B	2.0 mL	2.0 mL	2.0 mL	2.0 mL
Specimen	—	—	—	Serum 0.05 mL
Mix well and incubate at 37°C for 10 min.				
After cooling the reaction solution to room temperature, measure the absorbance of the test sample and standard solution with the blank solution as the control within the next 30min. Colorimeter with 550 nm filter Spectrophotometer : 550 nm				

*2 Sample blank is not required for common samples, but required for chyle or hemolyzed sample.

[Standard curve]

Microplate reader : SAFIRE (TECAN)

[Performance]

(1) Sensitivity

- The absorbance is below 0.07, when measuring purified water as a sample.
- The absorbance is 0.10 to 0.37, when measuring 1 mEq/L NEFA as a sample.

(2) Specificity

The NEFA concentration is less than ± 15% , when measuring the known concentration of control serum as a sample.

[Usage Notes]

(1) Sample

- Samples should be used immediately after collecting, because enzymes such as lipoprotein lipase, phospholipase *etc.* hydrolyze lipids and form fatty acids. Freeze sample, when a serum is stored.
- Hemolysis gives a higher value on the assay.
- Bilirubin gives a lower value on the assay.
- Ascorbic acid does not have significant effects on the assay.

(2) Notes on the assay

- Do not use the reagents described above in any procedures other than those described herein. Performance cannot be guaranteed if the reagents are used in other procedures or for other purposes.
- Operate the measurement apparatuses according to operator's manuals under appropriate conditions. Consult the apparatus manufacturer for details.
- Keep the provided reagents under the indicated conditions before the expiration.
- Do not use reagents, which were frozen in error. Such reagents may give false results.
- After opening the reagents, it is recommended to use them immediately. The opened vials should be capped and kept under the specified conditions.
- Do not use the containers and other materials in the package for any purpose other than those described herein.
- The vial is stoppered at reduced pressure. Slowly remove the stopper in order not to blow the powder in the vial.
- Do not use the reagents at other reaction temperatures and reac-

tion times than described herein.

- The reaction induced by Color Reagent A is completed in about 6 min at 37°C. Do not incubate more than 15 min.
- The reaction induced by Color Reagent B is completed in about 5 min at 37°C. When the incubation at 37°C is continued, absorbance decreases incrementally over time. Therefore, return it to room temperature immediately after the incubation for 10 min at 37°C.
- Avoid the direct sunlight during operation.
- This kit is for research use only. Not for diagnostic use.

(3) Safety precautions

- If a reagent comes into contact with the mouth, eyes, or skin, immediately wash with a lot of water. Consult a physician if necessary.
- A pipette with a safety pipette filler should be used.
- Be careful not to cut yourself with the aluminum cap when removing it from the vial.

(4) Waste

- The waste should be processed appropriately according to regulations.
- All the devices including reagents and vials that come into contact with the specimen should be considered potentially infectious.
- Chromogen Reagent A and Solvent A contain sodium azide as preservative (0.07% in Color A when reconstituted). NEFA Standard Solution contains 0.05% sodium azide. Sodium azide may react with lead or copper plumbing to form explosive compounds. Even though the reagent contains minute quantity of sodium azide, drains should be flushed well with a large amount of water, when discarding the reagents.

Expiration date : 12 months after the manufacture

Storage : Store at 2 ~ 10°C

Package : 750 tests

[Reference]

1. Shimizu, S., Yasui, K., Tani, Y. and Yamada, H. : *Biochem. Biophys. Res. Commun.*, **91**, 108 (1979).

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
U.S.A.
Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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