**FUJIFILM****Wako**

<For Research Use Only> Code No. 290-65901 (500 Tests)

## LabAssay™ Creatinine (Jaffé method)

### [Introduction]

Creatinine is a breakdown product of creatine phosphate in muscle. It is mainly filtered by the kidneys, though a small amount is actively secreted. LabAssay™ Creatinine is based on an *in vitro* colorimetric Jaffé method for the quantitative determination of creatinine in mouse serum or urine. It is a simultaneous multi-sample assay format using a microplate, but measurements can also be made using a test tube.

### [Kit contents]

(1)	Deproteinizing Reagent (Sodium Tungstate, Phosphoric Acid)	150mL	1vial
(2)	Picric Acid Reagent (Picric Acid 22mmol/L)	50mL	1vial
(3)	0.75mol/L Sodium Hydroxide Solution	50mL	1vial
(4)	Standard Solution (Creatinine 10mg/dL)	15mL	1vial

### [Materials and apparatuses to be prepared]

- 96wells microplate (transparent type)
- Micropipette
- Incubator maintaining at 25-30°C\*
- Plate mixer\*
- Microplate reader with 520nm wavelength filter  
(\* if the microplate reader is not equipped.)

### [For test tube method]

- Test tube
- Pipette
- Spectrophotometer or colorimeter with 500-530 nm wavelength filter

### [Assay principle]

The supernatant is separated, after adding deproteinizing reagent to a sample and centrifuging. Creatinine in the sample reacts with picric acid in alkaline solution to produce tangerine condensate (Jaffé reaction). Quantitation of total creatinine in the sample can be made by measurement of the absorbance.

### [Preparation of reagents to be used]

- ① Deproteinizing Reagent : ready-to-use
- ② Picric Acid Reagent : ready-to-use
- ③ 0.75mol/L Sodium Hydroxide Solution : ready-to-use

### ④ Standard solution

Standard solution is prepared by dilution of the provided Standard.

No.	Provided Std.	Distilled water for dilution	Sampling volume of diluted std.	Creatinine concentration
1	50 μL	150 μL	50 μL	2.5mg/dL
2	100 μL	100 μL	50 μL	5.0mg/dL
3	150 μL	50 μL	50 μL	7.5mg/dL
4	undiluted	—	50 μL	10.0mg/dL

### [Procedure]

#### (1) Assay in a microplate

Perform the assay according to the following table scheme.


	Test	Standard	Blank
	First operations should be performed in the sample tubes.		
Sample	Sample 50 μL	Standard solution 50 μL	Water 50 μL
Deproteinizing Reagent	300 μL	300 μL	300 μL
	Mix well, and stand at room temperature for 10 min. Centrifuge over 2,500rpm for 10 min. Subsequent operations should be performed in the wells.		
Supernatant	100 μL	100 μL	100 μL
	Stand at 25-30°C.		
Picric Acid Reagent	50 μL	50 μL	50 μL
0.75mol/L Sodium Hydroxide Solution	50 μL	50 μL	50 μL
	Stand at 25-30°C for 20 min. Measure the absorbance of the test sample and standard solution at 520nm with the blank solution as the control within the next 30 min.		

#### (2) Assay in a test tube

Perform the assay in a test tube according to the following table scheme.

	Test	Standard	Blank
Sample	Sample 0.5mL	Standard solution 0.5mL	Water 0.5mL
Deproteinizing Reagent	3.0mL	3.0mL	3.0mL
	Add with shaking, and stand at room temperature for 10 min. Centrifuge over 2,500rpm for 10 min.		
Supernatant	2.0mL	2.0mL	2.0mL
	Stand at 25-30°C bath.		
Picric Acid Reagent	1.0mL	1.0mL	1.0mL
0.75mol/L Sodium Hydroxide Solution	1.0mL	1.0mL	1.0mL
	Stand at 25-30°C for 20 min. Measure the absorbance of the test sample and standard solution with the blank solution as the control within the next 30 min. Colorimeter with 500-530nm wavelength filter Spectrophotometer : 520nm		

**[Standard curve]** [assay in a microplate]

Microplate reader : SAFIRE (TECAN)

**[Performance]**

(1) Sensitivity [assay in a test tube]

- The absorbance is 0.010 to 0.020, when measuring purified water as a sample.
- The absorbance is 0.400 to 0.500, when measuring 10mg/dL creatinine as a sample.

(2) Specificity

- The creatinine concentration is less than  $\pm 10\%$ , when measuring the known concentration of control serum as a sample.

**[Usage Notes]**

(1) Sample

- Heparin, citrate, oxalate and EDTA as an anticoagulant and sodium fluoride as an antiglycolysis may not affect the assay.
- Toluene should be used as the preservative for urine samples.
- A substance with an active methylene group causes an increase of the value.
- Ascorbic acid, bilirubin and hemolysis may not significantly affect the assay.

(2) Notes on the assay

- Do not use the reagents described above in any procedures other than those described herein. Performance cannot be guaranteed if the reagents are used in other procedures or for other purposes.
- Operate the measurement apparatuses according to operator's manuals under appropriate conditions. Consult the apparatus manufacturer for details.
- Keep the provided reagents under the indicated conditions before the expiration.
- The measurement should be done at an absorbance of 520nm, although maximum absorption of the pigment is 480nm, because picric acid may cause a high blank.
- Reaction temperature should be constant, because the temperature affects the assay value.
- Picric Acid Reagent and 0.75mol/L Sodium Hydroxide Solution can be mixed beforehand. However, the mixed solution should be used within the day.
- 0.75mol/L Sodium Hydroxide Solution should not be used if it is contaminated by Picric Acid Solution. Pay especially attentions not to contaminate with pipette.

- Urine creatinine samples below 10mg/dL should be measured without the dilution.
- Acetone should not be used for drying of tubes, because it effects the absorbance.
- In the case of urine, if there are no opacity after addition of Deproteinizing Reagent, it is not necessary to centrifuge.
- This kit is for research use only, but not for diagnostic use.

(3) Safety precautions

- If a reagent comes into contact with the mouth, eyes, or skin, immediately wash with a lot of water. Consult a physician if necessary.
- A pipette should be used with a safety pipette filler.
- Vials should be opened carefully.

(4) Waste

- The waste should be processed appropriately according to the law.
- All the devices including reagents and vials that come into contact with specimens should be considered potentially infectious.

**Expiration date** : 30 months after the manufacture

**Storage** : Store at 2 ~ 10°C

**Package** : 500 tests

**[References]**

1. Bonsnes, R. W. and Taussky, H. H. : *J. Biol. Chem.*, **158**, 581 (1945).
2. Henry, R. J. : *Clinical Chemistry*, 287 (Harper & Row), New York (1966).

---

**FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://www.wako-chem.co.jp>

**FUJIFILM Wako Chemicals U.S.A. Corporation**

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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